# ACOUTRAP

**PRODUCT SHEET** 





## **AUTOMATED HANDLING OF CELLS AND EXTRACELLULAR VESICLES**

The AcouTrap is a benchtop research instrument for sample preparation of cells, extracellular vesicles, bacteria and viruses from a range of different media.

The bioparticles are captured in an acoustic trap using gentle sound waves that keep the particles intact. The system is excellent for enrichment and purification as well as performing rapid reactions such as staining and antibody labeling. Mixing is enhanced by the acoustic streaming in the AcouTrap, thereby shortening incubation times.

The AcouTrap platform operates on standard 96-well microplates or microcentrifuge tubes. The system is automated and rapid, generating output samples with high quality, high throughput and high cell viability.

The trapping unit is available in two different designs, optimized for handling small to medium number of cells.

## **FEATURES**

- Automated sample processing
- Cell recovery > 90% with maintained viability
- Live viewing for process monitoring
- Compact benchtop design


Trapping unit	Small	Medium
Max. cell capacity*	100,000	500,000
Flow rate (µL/min)	15 – 100	50 – 300
Cell viability (%)**	99.3	98.9

<sup>\*</sup> Measured with cultured Jurkat cells. \*\* Compared to the input sample.





The central part of the system is an acoustic trapping unit where standing waves are used to trap and hold particles, cells or vesicles against a fluid flow. The acoustic forces can hold cells or polymer microbeads at sizes down to a few micrometers. For smaller particles, e.g. bacteria, extracellular vesicles or nanoparticles, larger polymer seed particles are initially trapped (A) and smaller objects are drawn to these particles by secondary acoustic forces (B). When the ultrasound is turned off (C) and the sample is ejected, and the particles dissociate (D).



## **TECHNICAL SPECIFICATIONS**

### **Dimensions**

 $51 \times 38 \times 37 \text{ cm}^3$  (width x depth x height)

## Weight

22 kg

## Power requirement

24 V DC, 5 A

## Fluid control

High precision syringe pump with 1 mL syringe

## External reservoirs for buffers, waste and cleaning agents

Five 125 mL bottles

### Microplate and microcentrifuge tube compatibility

Designed for use with standard 96-well microplates and 1.5 mL or 0.6 mL microcentrifuge tubes

#### **ACOUSORT**

AcouSort is an innovative medical technology company developing products for automated processing and refinement of biological samples such as blood or cell preparations.

## **Strategic Collaborations**

Our technology is available as OEM components, providing instrumentation manufacturers with a state-of-the-art ability to integrate sample processing steps that traditionally have had to be performed manually. AcouSort frequently engages in joint development projects with leading Life Science companies.

